

Module 2:

Basic MR Physics

Magnetic Resonance Imaging



An MR scanner consists of an electromagnet with a very strong magnetic field (1.5 - 7.0 Tesla)

Earth's magnetic field = 0.00005 Tesla

3 Tesla is 60,000 times stronger than the Earth's magnetic field.


What MRI Measures

- MRI is an extremely versatile imaging modality that can be used to study both brain structure and brain function.
- Both structural and functional MRI images are acquired using the same scanner.
- Different types of brain images can be generated to emphasize **contrast** related to different tissue characteristics.


Principles of MRI

- All magnetic resonance imaging techniques rely on a core set of physical principles.
- To understand we must begin by studying a **single atomic nuclei** and illustrate its impact on the generated MR signal.
- In particular we focus on **hydrogen** atoms consisting of a single proton.

Protons can be viewed as positively charged spheres which are always spinning. They give rise to a net magnetic moment along the axis of the spins.





Protons can be viewed as positively charged spheres which are always spinning. They give rise to a net magnetic moment along the axis of the spins.




Net Magnetization

- We cannot measure the magnetization of a single proton using MR, instead we measure the **net magnetization** of all nuclei within a volume.
- The net magnetization M can be viewed as a vector with two components.
 - A **longitudinal component** parallel to the magnetic field.
 - A **transverse component** perpendicular to the field.








In the absence of an external magnetic field, the nuclear magnetic moments are randomly oriented. There is no net magnetization.




When placed in a strong magnetic field, the nuclei align with the field. This creates a net longitudinal magnetization in the direction of the field.









The nuclei precess about the field with an angular frequency determined by the **Larmor frequency** but at a random phase.



A radio frequency (RF) pulse is used to align the phase and 'tip over' the nuclei. This causes the longitudinal magnetization to decrease, and establishes a new transversal magnetization.



A radio frequency (RF) pulse is used to align the phase and 'tip over' the nuclei. This causes the longitudinal magnetization to decrease, and establishes a new transversal magnetization.




Relaxation

- After the RF pulse is removed, the system seeks to return to equilibrium.
- The transverse magnetization starts to disappear (transversal relaxation), and the longitudinal magnetization grows back to its original size (longitudinal relaxation).
- During this process a signal is created that can be measured using a receiver coil.

Relaxation


- **Longitudinal Relaxation** is the restoration of net magnetization along the longitudinal direction as spins return to their parallel state.
 - Exponential growth described by time constant T_1
- **Transverse Relaxation** is the loss of net magnetization in the transverse plane due to loss of phase coherence.
 - Exponential decay described by time constant T_2

Longitudinal Relaxation Time




The restoration of longitudinal magnetization is described by a time constant T_1 .

Longitudinal Relaxation Time



Transverse Relaxation Time



The decay of magnetization due to interaction between nuclei is described by a time constant T_2 .

Image Contrast

- By altering how often we excite the nuclei (TR) and how soon after excitation we begin data collection (TE) we can control which characteristic is emphasized.
- The measured signal is approximately


$$M_0(1 - e^{-TR/T_1})e^{-TE/T_2}$$

where T_1 and T_2 are tissue properties.

Image Formation

- The goal of MRI is to construct an **image**, or a matrix of numbers that correspond to spatial locations.
- The image depicts the spatial distribution of some property of the nuclei within the sample.
- This could be the **density** of nuclei or the **relaxation time** of the tissues in which they reside.

Image Contrast



$$M_0(1 - e^{-TR/T_1})e^{-TE/T_2}$$

TE (echo time) -
the time between
excitation and
data collection.

fMRI Contrast

- T_2^* is the combined effect of T_2 and local inhomogeneities in the magnetic field.
- The scanner can be programmed to eliminate the effects of these inhomogeneities, or alternatively emphasize them.
- The latter types of procedures form the basis of BOLD fMRI.

Image Contrast

- Images can be produced that are sensitive primarily to T_1 , T_2 , or T_2^* .
- Because T_1 and T_2 vary with tissue type, they are able to represent boundaries between CSF, gray and white matter.
- Because T_2^* is sensitive to flow and oxygenation, it can be used to image brain function.

End of Module



@fMRIstats