

Module 3:

Image Formation

Signal Formation

- The subject is placed into the MR scanner.
 - Nuclei of ^1H atoms align with the magnetic field.
 - The nuclei precess about the field at similar frequencies, but at a random phase.
 - Net longitudinal magnetization in the direction of field.
- Within a slice, a radio frequency (RF) pulse is used to align the phase and ‘tip over’ the nuclei.
 - Causes the longitudinal magnetization to decrease, and establishes a new transversal magnetization.

Signal Formation

- After the RF pulse is removed, the system seeks to return to equilibrium.
 - The transverse magnetization disappears ([transversal relaxation](#)), and the longitudinal magnetization grows back to its original size ([longitudinal relaxation](#)).
 - Longitudinal relaxation: exponential growth described by time constant T_1 .
 - Transverse relaxation: exponential decay described by time constant T_2 .
- During this process a signal is created that can be measured using a receiver coil.

Slice Selection

- Most structural MRI and fMRI scans involve the construction of a three dimensional image from a set of two-dimensional slices.




Image Formation

- Imagine a brain slice split into a number of equally sized **volume elements** or **voxels**.





Image Formation

- Imagine a brain slice split into a number of equally sized **volume elements** or **voxels**.



Gradients

- The measured signal combines information from the whole brain:

$$S(t) = \iint \rho(x, y) dx dy$$

- A **magnetic field gradient** is used to sequentially control the spatial inhomogeneity of the magnetic field, so each measurement can be expressed:

$$S(k_x, k_y) = \iint \rho(x, y) e^{-i2\pi(k_x x + k_y y)} dx dy$$


K-space

- The measurements are acquired in the frequency-domain (**k-space**).
- By making measurements for multiple values of (k_x, k_y) we can gain enough information to solve the inverse problem and reconstruct $\rho(x, y)$.
- We can use the **inverse Fourier transform** (IFT):


$$\rho(x, y) = \iint S(k_x, k_y) e^{i2\pi(k_x x + k_y y)} dk_x dk_y$$

K-space Measurements

- In practice, data measurements are made **discretely** over a **finite region**.
 - Use discrete Fourier transforms.
- The number of k-space measurements we make influences the spatial resolution of the image.
 - Need enough measurements to solve inverse problem.





16 unknowns



4 unknowns

EPI and Spirals



Magnitude Images



- The measured k-space data is complex valued.
- We typically work with magnitude images, or

$$|\rho(x)| = \sqrt{\rho_R(x)^2 + \rho_I(x)^2}$$

where ρ_R and ρ_I are the real and imaginary parts of the k-space measurement.



Image Formation

k-space



IFT
→
←
FT

Image space



End of Module



@fMRIstats